


Anti-Cd5 Monoclonal Antibody () For Flow Cytometry - P06127

Item Number: CM0007989



Introduction

Rabbit anti-Human CD5 monoclonal antibody validated for flow cytometry (FC). Targets T-cell surface glycoprotein CD5 (T1/Leu-1), a key immune regulator. Detects human CD5, predicted molecular weight ~55 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD5 mAb
Remarks/Alias	T1; LEU1
Species	Human
GeneID (Human)	921
GeneID	921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of human CD5 (NP_055022.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag