


Anti-Cd5 Rabbit Monoclonal Antibody For Fc, Elisa - P06127

Item Number: CM0009647



Introduction

Rabbit monoclonal antibody targeting human CD5 (T1/Leu-1), validated for flow cytometry (FC) and ELISA. Ideal for studying T-cell signaling, thymocyte selection, and B1a cell biology. Research-use only.

[Learn More](#)

Parameter	Value
Product Name	CD5 Rabbit mAb
Remarks/Alias	T1; LEU1
Species	Human
Gene ID (Human)	921
Gene ID	921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of human CD5 (NP_055022.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag