


Anti-Cd5 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06127

Item Number: CM0011235



Introduction

High-quality CD5 Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human CD5, a T-cell and B1a cell surface receptor involved in immune regulation. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD5 Rabbit mAb
Remarks/Alias	T1; LEU1; CD5
Species	Human
GeneID (Human)	921
GeneID	921
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD5 (P06127).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag