


Anti-Cd5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P06127

Item Number: CM0030951



Introduction

CD5 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Demonstrates cross-reactivity with human, mouse, and rat CD5. Ideal for T-cell signaling and B-cell interaction studies.

[Learn More](#)

Parameter	Value
Product Name	CD5 Rabbit pAb
Remarks/Alias	T1; LEU1; CD5
Species	Human
GeneID (human)	921
GeneID	921
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 158-371 of human CD5 (NP_055022.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag