

Anti-Cd50 (Icam-3) Monoclonal Antibody For Flow Cytometry - P32942

Item Number: CM0009981



Introduction

Rabbit monoclonal antibody against human CD50 (ICAM-3), clone . Validated for flow cytometry (FC). Immunogen: recombinant protein (aa 30-114). Suitable for leukocyte adhesion and immune interaction studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD50(ICAM-3) mAb
Remarks/Alias	CD50; CDW50; ICAM-R
Species	Human
GeneID (human)	3385
GeneID	3385
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-114 of human CD50(ICAM-3) (NP_002153.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32942
Protein Weight	60kDa
Shipping	Ice bag