

Anti-Cd52 Polyclonal Antibody For Wb - P31358

Item Number: CM0022090



Introduction

Rabbit polyclonal antibody specific to human CD52 (CAMPATH-1 antigen), a 7 kDa GPI-anchored membrane glycoprotein. Validated for WB, IHC-P, ELISA, and cross-reacts with mouse. Ideal for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CD52 Rabbit pAb
Remarks/Alias	HE5; CDW52; EDDM5; CD52
Species	Human
GeneID (human)	1043
GeneID	1043
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-61 of human CD52 (NP_001794.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P31358
Protein Weight	7kDa
Shipping	Ice bag