

Anti-Cd53 Rabbit Polyclonal Antibody For Wb, Elisa - P19397

Item Number: CM0018170



Introduction

High-quality rabbit polyclonal antibody targeting human CD53 (tetraspanin-25), validated for WB and ELISA. Cross-reacts with mouse. Ideal for immune cell signaling and B/T cell research.

[Learn More](#)

Parameter	Specification
Product Name	CD53 Rabbit pAb
Remarks/Alias	MOX44; TSPAN25; CD53
Species	Human
Gene ID (Human)	963
Gene ID	963
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human CD53 (NP_000551.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19397
Protein Weight	24kDa
Shipping	Ice bag