


Anti-Cd55 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P08174

Item Number: CM0003220



Introduction

Anti-CD55 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA in Human. Detects complement decay-accelerating factor (DAF, CD55) at 41kDa. Ideal for complement regulation and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	CD55 Rabbit mAb
Remarks/Alias	CR; TC; DAF; CROM; CHAPLE; CD55
Species	Human
GeneID (Human)	1604
GeneID	1604
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human CD55 (P08174).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08174
Protein Weight	41kDa
Shipping	Ice bag