


Anti-Cd55 Rabbit Polyclonal Antibody For Wb, Elisa - P08174

Item Number: CM0015738



Introduction

High-quality rabbit polyclonal antibody against human CD55 (DAF), validated for Western blot and ELISA. Suitable for human and mouse samples. Ideal for complement research and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	CD55 Rabbit pAb
Remarks/Alias	CR; TC; DAF; CROM; CHAPLE; CD55
Species	Human
GeneID (Human)	1604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human CD55 (NP_000565.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08174
Protein Weight	41kDa
Shipping	Ice bag