


Anti-Cd55 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08174

Item Number: CM0014520



Introduction

CD55 Rabbit Polyclonal Antibody (DAF) validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets complement decay-accelerating factor, a key regulator of innate immunity and viral entry receptor.

[Learn More](#)

Parameter	Value
Product Name	CD55 Rabbit pAb
Remarks/Alias	CR, TC, DAF, CROM, CHAPLE, CD55
Species	Human
Gene ID (Human)	1604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-300 of human CD55 (NP_000565.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08174
Protein Weight	41kDa
Shipping	Ice bag