


Anti-Cd55 Rabbit Polymab® For Wb, Elisa - P08174

Item Number: CM0007440



Introduction

CD55 (DAF) Rabbit PolymAb® validated for Western blot and ELISA; cross-reacts with human and mouse. Recombinant monoclonal antibody targeting complement decay-accelerating factor, a key regulator of complement cascade and viral entry receptor.

[Learn More](#)

Parameter	Value
Product Name	CD55 Rabbit PolymAb®
Remarks/Alias	CR; TC; DAF; CROM; CHAPLE
Species	Mouse/Human
GeneID (Human)	1604
GeneID	13136/1604
Immunogen	Recombinant protein of human CD55
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08174
Protein Weight	41kDa
Shipping	Ice bag