

Anti-Cd56/Ncam1 Monoclonal Antibody For Fc - P13591

Item Number: CM0012584



Introduction

Recombinant rabbit monoclonal antibody against human CD56/NCAM1, validated for flow cytometry (FC). Targets the neural cell adhesion molecule involved in neuron adhesion and viral entry receptor (rabies, Zika).

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Human CD56/NCAM1 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
GeneID (Human)	4684
GeneID	4684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-413 of human CD56/NCAM1 (NP_851996.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13591
Protein Weight	95kDa
Shipping	Ice bag