


Anti-Cd58 (Lfa-3) Rabbit Monoclonal Antibody For Wb, Elisa - P19256

Item Number: CM0006559



Introduction

Recombinant rabbit monoclonal antibody against human CD58 (LFA-3), validated for Western blot and ELISA. Recognizes endogenous CD58 at ~28 kDa. Suitable for studying T-cell adhesion and immunological synapse formation.

[Learn More](#)

Parameter	Value
Product Name	CD58 Rabbit mAb
Remarks/Alias	ag3; LFA3; LFA-3; CD58
Species	Human
Gene ID (Human)	965
Gene ID	965
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-215 of human CD58 (NP_001770.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P19256
Protein Weight	28 kDa
Shipping	Ice bag