


Anti-Cd59 Polyclonal Antibody For Wb, Elisa - P13987

Item Number: CM0013793



Introduction

Rabbit polyclonal antibody specific for human CD59, validated for WB and ELISA. CD59 is a GPI-anchored glycoprotein that inhibits complement MAC formation, a key target in hemolytic anemia research.

[Learn More](#)

Parameter	Value
Product Name	CD59 Rabbit pAb
Remarks/Alias	1F5; EJ16; EJ30; EL32; G344; MIN1; MIN2; MIN3; M1RL; HRF20; MAC1F; MEM43; MIC11; MSK21; 16.3A5; HRF-20; MAC-IP; p18-20; CD59
Species	Human
Gene ID (Human)	966
Gene ID	966
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-128 of human CD59 (NP_001120695.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13987
Protein Weight	14kDa
Shipping	Ice bag