


Anti-Cd59 Rabbit Monoclonal Antibody For Wb, Elisa - P13987

Item Number: CM0010989



Introduction

Recombinant rabbit monoclonal antibody against human CD59 (protectin) for WB and ELISA. Suitable for complement regulation studies and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD59 Rabbit mAb
Remarks/Alias	1F5; EJ16; EJ30; EL32; G344; MIN1; MIN2; MIN3; MIRL; HRF20; MACIF; MEM43; MIC11; MSK21; 16.3A5; HRF-20; MAC-IP; p18-20; CD59
Species	Human
GeneID (Human)	966
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 49-128 of human CD59 (NP_000602.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13987
Protein Weight	14kDa
Shipping	Ice bag