


Anti-Cd59 Rabbit Monoclonal Antibody For Wb, Elisa - P13987

Item Number: CM0006851



Introduction

Rabbit monoclonal antibody against human CD59 validated for Western blot and ELISA. Targets the complement membrane attack complex inhibitor CD59 glycoprotein, suitable for complement research and diagnostic assays.

[Learn More](#)

Parameter	Specification
Product Name	CD59 Rabbit mAb
Remarks/Alias	1F5; EJ16; EJ30; EL32; G344; MIN1; MIN2; MIN3; M1RL; HRF20; MACIF; MEM43; MIC11; MSK21; 16.3A5; HRF-20; MAC-IP; p18-20; CD59
Species	Human
GeneID (human)	966
GeneID	966
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 26-126 of human CD59. (NP_000602.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13987
Protein Weight	14 kDa
Shipping	Ice bag