


Anti-Cd6 Rabbit Polyclonal Antibody For Wb & Elisa - P30203

Item Number: CM0028918



Introduction

High-quality rabbit polyclonal antibody targeting human CD6 (TP120). Suitable for Western blot and ELISA applications. Ideal for T-cell biology, immunology, and infection research. Cross-reacts with human CD6.

[Learn More](#)

Parameter	Value
Product Name	CD6 Rabbit pAb
Remarks/Alias	TP120; CD6
Species	Human
Gene ID (Human)	923
Gene ID	923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 424-668 of human CD6 (NP_006716.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P30203
Protein Weight	72kDa
Shipping	Ice bag