

Anti-Cd62E Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - Q00690

Item Number: CM0011482

Introduction

Rabbit monoclonal antibody against mouse CD62E (E-selectin), validated for WB, IF/ICC, FC, and ELISA. Recombinant immunogen targeting residues 22-557. Reacts with mouse, 66 kDa. Ideal for endothelial adhesion and leukocyte trafficking research.

[Learn More](#)

Parameter	Value
Product Name	CD62E Rabbit mAb
Remarks/Alias	Elam; CD62E; ELAM-1; LECAM2; E-selectin
Species	Mouse
Gene ID (Human)	6401
Gene ID	20339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-557 of mouse CD62E (NP_035475.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Mouse
SWISS	Q00690
Protein Weight	66kDa
Shipping	Ice bag