

Anti-Cd62L/L-Selectin Monoclonal Antibody (Apc/Cyanine7) For Flow Cytometry - P14151

Item Number: CM0009452



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody against human CD62L (L-selectin). Validated for flow cytometry. Ideal for studying lymphocyte homing, leukocyte rolling, and cell adhesion in human samples.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human CD62L/L-Selectin mAb
Remarks/Alias	TQ1; LAM1; LEU8; LNHR; LSEL; CD62L; LYAM1; PLNHR; LECAM1
Species	Human
GeneID	6402
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of human CD62L/L-Selectin (NP_000646.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P14151
Protein Weight	42kDa
Shipping	Ice bag