

Anti-Cd62L/L-Selectin Polyclonal Antibody For Wb, Ihc-P, Elisa - P14151

Item Number: CM0023516



Introduction

Rabbit polyclonal antibody against human CD62L/L-Selectin, suitable for WB, IHC-P, and ELISA. Reacts with human, mouse, rat. Ideal for leukocyte adhesion and lymphocyte homing research.

[Learn More](#)

Parameter	Value
Product Name	CD62L/L-Selectin Rabbit pAb
Remarks/Alias	TQ1; LAM1; LEU8; LNHR; LSEL; CD62L; LYAM1; PLNHR; LECAM1; CD62L/L-Selectin
Species	Human
GeneID (Human)	6402
GeneID	6402
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD62L/L-Selectin (NP_000646.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14151
Protein Weight	42kDa
Shipping	Ice bag