

## Anti-Cd62P/P-Selectin Polyclonal Antibody For Wb, Elisa - P16109

Item Number: CM0015517



### Introduction

CamelBio's rabbit anti-human CD62P/P-selectin polyclonal antibody targets the extracellular domain (aa 220-319). Suitable for WB and ELISA. Cross-reacts with human and mouse P-selectin (91kDa). Supports research on inflammation, thrombosis, and cell adhesion.

[Learn More](#)

| Parameter        | Value                                                                                                                         |
|------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD62P/P-selectin Rabbit pAb                                                                                                   |
| Remarks/Alias    | CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM; CD62P/P-selectin                                                             |
| Species          | Human                                                                                                                         |
| GeneID (Human)   | 6403                                                                                                                          |
| GeneID           | 6403                                                                                                                          |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 220-319 of human CD62P/P-selectin (NP_002996.2) |
| Source           | Rabbit                                                                                                                        |
| Category         | Polyclonal Antibodies                                                                                                         |
| Application      | WB, ELISA                                                                                                                     |
| Cross Reactivity | Human, Mouse                                                                                                                  |
| SWISS            | P16109                                                                                                                        |
| Protein Weight   | 91kDa                                                                                                                         |
| Shipping         | Ice bag                                                                                                                       |