

## Anti-Cd62P/P-Selectin Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P16109

Item Number: CM0007041



### Introduction

This CD62P/P-selectin Rabbit Monoclonal Antibody (mAb) is validated for WB, IF/ICC, and ELISA in human samples. It targets P-selectin, a key adhesion molecule involved in leukocyte rolling and inflammation.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD62P/P-selectin Rabbit mAb                                                                                       |
| Remarks/Alias    | CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM; CD62P/P-selectin                                                 |
| Species          | Human                                                                                                             |
| Gene ID          | 6403 (Human)                                                                                                      |
| Immunogen        | Recombinant Protein corresponding to a sequence within amino acids 42-771 of human CD62P/P-selectin (NP_002996.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, ELISA                                                                                                 |
| Cross Reactivity | Human                                                                                                             |
| SWISS            | P16109                                                                                                            |
| Protein Weight   | 91kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |