

Anti-Cd62P/P-Selectin Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P16109

Item Number: CM0007796



Introduction

Recombinant rabbit monoclonal antibody against human CD62P/P-selectin, validated for WB, IHC-P, IF-P, and ELISA. Ideal for studying leukocyte adhesion and inflammation signaling pathways in research use.

[Learn More](#)

Parameter	Value
Product Name	CD62P/P-selectin Rabbit mAb
Remarks/Alias	CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM
Species	Human
GeneID (Human)	6403
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P16109
Protein Weight	91kDa
Shipping	Ice bag