

Anti-Cd62P/P-Selectin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P16109

Item Number: CM0018486



Introduction

CD62P/P-selectin rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human and rat. Ideal for inflammation, leukocyte adhesion, and platelet activation research from CamelBio's IVD raw material platform.

[Learn More](#)

Parameter	Value
Product Name	CD62P/P-selectin Rabbit pAb
Remarks/Alias	CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM
Species	Human
Gene ID (Human)	6403
Gene ID	6403
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-200 of human CD62P/P-selectin (NP_002996.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P16109
Protein Weight	91kDa
Shipping	Ice bag