


Anti-Cd63 Monoclonal Antibody For Wb, Ihc-P, Elisa - P08962

Item Number: CM0003142



Introduction

CD63 Rabbit mAb validated for WB, IHC-P, ELISA. Targets human CD63 (granulophysin, LAMP-3), a tetraspanin involved in cell adhesion, signaling, and exosome biology. Ideal for cancer and exosome research.

[Learn More](#)

Parameter	Value
Product Name	CD63 Rabbit mAb
Remarks/Alias	MLA1, ME491, LAMP-3, OMA81H, TSPAN30, CD63
Species	Human
GeneID (Human)	967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 103-203 of human CD63 (NP_001771.1)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08962
Protein Weight	17kDa/23kDa/25kDa
Shipping	Ice bag