


Anti-Cd63 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08962

Item Number: CM0026817



Introduction

High-quality rabbit polyclonal antibody targeting human CD63, suitable for Western blot, IHC-P, IF/ICC, and ELISA. Ideal for studying exosome biogenesis, melanoma, and tetraspanin functions. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CD63 Rabbit pAb
Remarks/Alias	MLA1; ME491; LAMP-3; OMA81H; TSPAN30; CD63
Species	Human
GeneID (Human)	967
GeneID	967
Immunogen	Recombinant protein corresponding to a sequence containing amino acids 103-203 of human CD63 (NP_001771.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P08962
Protein Weight	26kDa
Shipping	Ice bag