

Anti-Cd64 Rabbit Monoclonal Antibody For Wb, Elisa - P12314

Item Number: CM0006679



Introduction

Rabbit monoclonal antibody targeting human CD64 (FcγRI), validated for Western blot and ELISA applications. High-affinity IgG Fc receptor selectively expressed on monocytes/macrophages, mediating ADCC.

[Learn More](#)

Parameter	Value
Product Name	CD64 Rabbit mAb
Remarks/Alias	CD64; FCG1; FCRI; CD64A; FCGR1; IGFR1; FcgammaRI
Species	Human
Gene ID (Human)	2209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-288 of human CD64 (NP_000557.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P12314
Protein Weight	45-75kDa
Shipping	Ice bag