

Anti-Cd64 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P12314

Item Number: CM0014078



Introduction

High-affinity IgG Fc receptor rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects CD64 (FCGR1A) in human, mouse, and rat. Suitable for immunology, ADCC, and diagnostic immunoassay development.

[Learn More](#)

Parameter	Value
Product Name	CD64 Rabbit pAb
Remarks/Alias	CD64; FCG1; FCRI; CD64A; FCGR1; IGFR1; FcgammaRI
Species	Human
GeneID (Human)	2209
GeneID	2209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-288 of human CD64 (NP_000557.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12314
Protein Weight	43kDa
Shipping	Ice bag