

Anti-Cd66A/Ceacam1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Fc, Elisa - P31809

Item Number: CM0009489



Introduction

Rabbit monoclonal antibody targeting mouse CD66a/CEACAM1, validated for WB, IF-P, IHC-P, FC, and ELISA, with rat cross-reactivity. Ideal for investigating immune regulation, cell adhesion, and coronavirus infection research.

[Learn More](#)

Parameter	Value
Product Name	CD66a/CEACAM1 Rabbit mAb
Remarks/Alias	Bgp; Cc1; Hv2; Bgp1; Cea1; Cea7; Hv-2; MHVR; bb-1; C-CAM; CD66a; Cea-1; Cea-7; MHVR1; Mhv-1; mCEA1; mmCGM1; mmCGM2; mmCGM1a
Species	Mouse
GeneID (Human)	634
GeneID	26365
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-428 of mouse CD66a/CEACAM1a (NP_001034274.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, FC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P31809
Protein Weight	57kDa/50kDa/37kDa/30kDa
Shipping	Ice bag