

Anti-Cd66A/Ceacam1 Rabbit Polymab® For Wb, Ihc, If, Elisa - P13688 / P31809

Item Number: CM0012791



Introduction

Rabbit PolymAb® targeting human/mouse CD66a/CEACAM1, validated for WB, IHC, IF, ELISA. Mediates cell adhesion and immune regulation; ideal for cancer and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD66a/CEACAM1 Rabbit PolymAb®
Remarks/Alias	Bgp; Cc1; Hv2; Bgp1; Cea1; Cea7; Hv-2; MHVR; bb-1; C-CAM; CD66a; Cea-1; Cea-7; MHVR1; Mhv-1; mCEA1; mmCGM1; mmCGM2; mmCGM1a
Species	Human/Mouse
GeneID (Human)	634
GeneID	634/26365
Immunogen	Recombinant protein of human/mouse CD66a/CEACAM1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13688/P31809
Protein Weight	30kDa/37kDa/43kDa/46kDa/50kDa/57kDa
Shipping	Ice bag