


## Anti-Cd66B/Ceacam8 Polyclonal Antibody For Wb, If/Icc, Elisa - P31997

Item Number: CM0026483



### Introduction

Rabbit polyclonal antibody targeting human CD66b/CEACAM8 (38 kDa) for Western blot, IF/ICC, and ELISA. Suitable for investigating neutrophil cell adhesion and chronic myeloid leukemia research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD66b/CEACAM8 Rabbit pAb                                                                                                  |
| Remarks/Alias    | CD67; CGM6; CD66b; NCA-95                                                                                                 |
| Species          | Human                                                                                                                     |
| GeneID (Human)   | 1088                                                                                                                      |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 35-320 of human CD66b/CEACAM8 (NP_001807.2) |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IF/ICC, ELISA                                                                                                         |
| Cross Reactivity | Human                                                                                                                     |
| SWISS            | P31997                                                                                                                    |
| Protein Weight   | 38kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |