

Anti-Cd66D/Ceacam3 Polyclonal Antibody For Wb, Elisa - P40198

Item Number: CM0027836



Introduction

Rabbit polyclonal antibody targeting CD66d/CEACAM3, a granulocyte receptor mediating phagocytosis of pathogens. Validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	CD66d/CEACAM3 Rabbit pAb
Remarks/Alias	CEA; CGM1; W264; W282; CD66D; CD66d/CEACAM3
Species	Human
Gene ID (Human)	1084
Gene ID	1084
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human CD66d/CEACAM3 (NP_001806.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P40198
Protein Weight	27kDa
Shipping	Ice bag