

Anti-Cd68 Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - P31996

Item Number: CM0006627

Introduction

Rabbit monoclonal antibody targeting mouse CD68 (Macrosialin). Validated for WB, IF/ICC, FC, ELISA. Recognizes ~33/34 kDa protein. Suitable for studying phagocytosis, macrophage homing, and lysosomal metabolism.

[Learn More](#)

Parameter	Value
Product Name	CD68 Rabbit mAb
Remarks/Alias	Lamp4; gp110; Scard1; CD68
Species	Mouse
GeneID (Human)	968
GeneID (Mouse)	12514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-291 of mouse CD68 (NP_001277987.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Mouse
SWISS	P31996
Protein Weight	33KDa/34KDa
Shipping	Ice bag