

Anti-Cd68 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P34810

Item Number: CM0004185

Introduction

High-quality anti-CD68 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Recognizes human macrosialin (37kDa), a macrophage marker involved in phagocytosis and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	CD68 Rabbit mAb
Remarks/Alias	GP110; LAMP4; SCARD1; CD68
Species	Human
GeneID (Human)	968
GeneID	968
Immunogen	Synthetic peptide corresponding to amino acids 100-200 of human CD68 (NP_001242.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P34810
Protein Weight	37kDa
Shipping	Ice bag