

Anti-Cd68 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P34810

Item Number: CM0025441

Introduction

High-quality rabbit polyclonal antibody against human CD68 (GP110, LAMP4, SCARD1) for western blot, immunofluorescence, and ELISA. Validated for human, mouse, and rat reactivity. Suitable for macrophage and tumor studies.

[Learn More](#)

Parameter	Value
Product Name	CD68 Rabbit pAb
Remarks/Alias	GP110; LAMP4; SCARD1; CD68
Species	Human
GeneID (Human)	968
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD68 (NP_001242.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34810
Protein Weight	37kDa
Shipping	Ice bag