

Anti-Cd68 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P31996

Item Number: CM0022959

Introduction

CD68 rabbit polyclonal antibody (Lamp4/gp110/Scard1) validated for WB, IF/ICC, IHC-P, ELISA in human, mouse, rat. Ideal for macrophage and phagocytosis studies.

[Learn More](#)

Parameter	Value
Product Name	CD68 Rabbit pAb
Remarks/Alias	Lamp4; gp110; Scard1; CD68
Species	Mouse
GeneID (Human)	968
GeneID	12514
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of mouse CD68 (NP_001277987.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31996
Protein Weight	35 kDa
Shipping	Ice bag