

## Anti-Cd68 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P34810

Item Number: CM0013967

### Introduction

Rabbit polyclonal antibody against human CD68 (macrosialin, Gp110) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse. Suitable for macrophage and phagocytosis research.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | CD68 Rabbit pAb                                                                                     |
| Remarks/Alias    | GP110; LAMP4; SCARD1; CD68                                                                          |
| Species          | Human                                                                                               |
| Gene ID          | 968                                                                                                 |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 22-319 of human CD68 (NP_001242.2) |
| Source           | Rabbit                                                                                              |
| Category         | Polyclonal Antibodies                                                                               |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                            |
| Cross Reactivity | Human, Mouse                                                                                        |
| SWISS            | P34810                                                                                              |
| Protein Weight   | 37kDa                                                                                               |
| Shipping         | Ice bag                                                                                             |