

Anti-Cd68 Rabbit Polymab® For Wb, Ihc-P, If/Icc, Elisa - P34810 / P31996

Item Number: CM0005427

Introduction

Recombinant rabbit monoclonal CD68 antibody validated for WB, IHC-P, IF/ICC, ELISA.

Recognizes human and mouse macrosialin (Lamp4, gp110, Scard1), a

lysosomal/endosomal membrane glycoprotein involved in macrophage phagocytic activities and tumor cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	CD68 Rabbit PolymAb®
Remarks/Alias	Lamp4; gp110; Scard1
Species	Human/Mouse
GeneID (Human)	968
GeneID	968/12514
Immunogen	Recombinant protein of human/mouse CD68
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P34810/P31996
Protein Weight	70-130 kDa
Shipping	Ice bag