

Anti-Cd69 Monoclonal Antibody For Fc - Q07108

Item Number: CM0007789



Introduction

Rabbit anti-Human/Monkey CD69 mAb for flow cytometry. Recombinant protein immunogen, cross-reacts with human and cynomolgus. Ideal for early T-cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD69 mAb
Remarks/Alias	AIM, EA1, MLR-3, CLEC2C, GP32/28, BL-AC/P26
Species	Human
Gene ID (Human)	969
Gene ID	969
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 78-199 of human CD69 (NP_001772.1).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q07108
Protein Weight	23kDa
Shipping	Ice bag