

Anti-Cd69 Rabbit Polymab® For Wb, Ihc-P, If-P, Elisa - Q07108

Item Number: CM0012798



Introduction

Rabbit monoclonal antibody against human CD69 (AIM, EA1, CLEC2C) for WB, IHC-P, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for T-cell activation and immune homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	CD69 Rabbit PolymAb®
Remarks/Alias	AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26
Species	Human
GeneID	969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 62-199 of human CD69 (NP_001772.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07108
Protein Weight	23kDa
Shipping	Ice bag