

Anti-Cd7 Rabbit Monoclonal Antibody For Wb, Elisa - P09564

Item Number: CM0007611

Introduction

High-quality CD7 Rabbit mAb for WB and ELISA applications. Specific for human CD7, a T-cell and NK cell marker involved in immune activation. Ideal for immunology research and diagnostic assay development.

[Learn More](#)

Parameter	Specification
Product Name	CD7 Rabbit mAb
Remarks/Alias	GP40; TP41; LEU-9; CD7
Species	Human
Gene ID (Human)	924
Gene ID	924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7. (NP_006128.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag