

Anti-Cd7 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P09564

Item Number: CM0008813

Introduction

High-quality rabbit monoclonal antibody against human CD7, validated for WB, IHC-P, and ELISA. CD7 is a T-cell and NK cell transmembrane glycoprotein involved in adaptive immunity; ideal for immunology research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	CD7 Rabbit mAb
Remarks/Alias	GP40; TP41; Tp40; LEU-9; CD7
Species	Human
Gene ID (Human)	924
Gene ID	924
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 141-240 of human CD7 (P09564).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag