


Anti-Cd70 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P32970

Item Number: CM0019347



Introduction

High-quality rabbit polyclonal antibody against human CD70 (CD27 ligand, TNFSF7) validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for immunology research on T cell activation and B cell immune surveillance.

[Learn More](#)

Parameter	Value
Product Name	CD70 Rabbit pAb
Remarks/Alias	CD27L; LPFS3; CD27-L; CD27LG; TNFSF7; TNLG8A; CD70
Species	Human
GeneID (Human)	970
GeneID	970
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-126 of human CD70 (NP_001243.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32970
Protein Weight	21kDa
Shipping	Ice bag