


Anti-Cd71 Monoclonal Antibody () For Flow Cytometry - P02786

Item Number: CM0009424



Introduction

-conjugated rabbit monoclonal antibody targeting human CD71 (transferrin receptor 1). Validated for flow cytometry, ideal for studying iron metabolism, T/B cell proliferation, and viral entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD71 mAb
Remarks/Alias	T9; TR; TFR; p90; CD71; TFR1; TRFR; IMD46
Species	Human
GeneID (Human)	7037
GeneID	7037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-250 of human CD71 (NP_003225.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P02786
Protein Weight	38kDa
Shipping	Ice bag