


Anti-Cd71/Transferrin Receptor Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02786

Item Number: CM0006695



Introduction

KO-validated rabbit monoclonal antibody against human CD71/Transferrin Receptor (TfR1), suitable for WB, IHC-P, IF/ICC, and ELISA. Ideal for iron metabolism, oncology, and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CD71/Transferrin Receptor Rabbit mAb
Remarks/Alias	T9; TR; TFR; p90; CD71; TFR1; TRFR; IMD46
Species	Human
GeneID (Human)	7037
GeneID	7037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-250 of human CD71/Transferrin Receptor (NP_003225.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P02786
Protein Weight	85kDa
Shipping	Ice bag