


Anti-Cd71/Transferrin Receptor Rabbit Monoclonal Antibody For Wb - P02786

Item Number: CM0005673



Introduction

Rabbit monoclonal antibody against human CD71/Transferrin Receptor (TFRC). Validated for WB, IP, IF/ICC, IF-P, IHC-P, ELISA. Reacts with human and mouse. Ideal for iron metabolism, viral entry, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CD71/Transferrin Receptor Rabbit mAb
Remarks/Alias	T9, TR, TFR, p90, CD71, TFR1, TRFR, IMD46
Species	Human
GeneID (Human)	7037
GeneID	7037
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02786
Protein Weight	85 kDa
Shipping	Ice bag