


Anti-Cd71/Transferrin Receptor Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02786

Item Number: CM0027018



Introduction

Rabbit polyclonal antibody against human CD71/Transferrin Receptor (TFRC, TFR1, p90). Validated for WB, IHC-P, IF/ICC, and ELISA. Relevant for iron metabolism and viral entry research - P02786.

[Learn More](#)

Parameter	Value
Product Name	CD71/Transferrin Receptor Rabbit pAb
Remarks/Alias	T9; TR; TFR; p90; CD71; TFR1; TRFR; IMD46; CD71/Transferrin Receptor
Species	Human
GeneID	7037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human CD71/Transferrin Receptor (NP_001121620.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P02786
Protein Weight	85kDa
Shipping	Ice bag