

Anti-Cd72 Rabbit Polyclonal Antibody For Wb, Elisa - P21854

Item Number: CM0025586

Introduction

High-quality rabbit anti-CD72 polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat CD72, a B-cell co-receptor regulating BCR signaling and immune homeostasis. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD72 Rabbit pAb
Remarks/Alias	LYB2, CD72b, CD72
Species	Human
GeneID (Human)	971
GeneID	971
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 120-359 of human CD72 (NP_001773.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21854
Protein Weight	40kDa
Shipping	Ice bag