

Anti-Cd73/Nt5E Polyclonal Antibody For Wb, Ihc-P, And Elisa - P21589

Item Number: CM0020884



Introduction

Anti-CD73/NT5E rabbit polyclonal antibody validated for Western blot, IHC paraffin, and ELISA. Recognizes human, mouse, and rat CD73 ecto-5'-nucleotidase. Suitable for immunodetection of this GPI-anchored cell membrane enzyme involved in nucleotide hydrolysis.

[Learn More](#)

Parameter	Value
Product Name	CD73/NT5E Rabbit pAb
Remarks/Alias	NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA; CD73/NT5E
Species	Human
GeneID (Human)	4907
GeneID	4907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-264 of human CD73/NT5E (NP_002517.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB,IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21589
Protein Weight	63kDa
Shipping	Ice bag