

Anti-Cd73/Nt5E Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P21589

Item Number: CM0020390



Introduction

Rabbit polyclonal antibody targeting human CD73/NT5E (ecto-5'-nucleotidase), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suits studies in nucleotide metabolism and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	CD73/NT5E Rabbit pAb
Remarks/Alias	NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA; 5E
Species	Human
GeneID (Human)	4907
GeneID	4907
Immunogen	Recombinant protein (recombinant fusion protein containing amino acids 27-377 of human CD73/NT5E, NP_002517.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21589
Protein Weight	63kDa
Shipping	Ice bag